


MINT1 Antibody / APBA1 (R33717)

Catalog No.	Formulation	Size
R33717-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	320
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This MINT1 antibody is available for research use only.



Western blot testing of human cerebellum lysate with MINT1 antibody at 0.1ug/ml.
Predicted molecular weight: ~93kDa.

Description

Additional name(s) for this target protein: Amyloid beta (A4) precursor protein-binding, family A, member 1; APBA1

Application Notes

Optimal dilution of the MINT1 antibody should be determined by the researcher.

Immunogen

Amino acids EASHPSQDGKRQYK were used as the immunogen for this MINT1 antibody.

Storage

Aliquot and store the MINT1 antibody at -20oC.