


MIPEP Antibody / Mitochondrial intermediate peptidase [clone 166CT22.3.7] (F40233)

Catalog No.	Formulation	Size
F40233-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40233-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, k
Clone Name	166CT22.3.7
Purity	Purified
UniProt	Q99797
Applications	Western Blot : 1:100-1:1000
Limitations	This MIPEP antibody is available for research use only.



Description

Cleaves proteins, imported into the mitochondrion, to their mature size.

Application Notes

Titration of the MIPEP antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This MIPEP antibody was produced from a mouse immunized with MIPEP recombinant protein.

Storage

Aliquot the MIPEP antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.