


Mkl Antibody (F43685)


Catalog No.	Formulation	Size
F43685-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F43685-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9D2Y4
Applications	Western Blot : 1:1000
Limitations	This Mkl antibody is available for research use only.



Western blot testing of Mkl antibody at 1:2000 dilution and mouse samples: (1) lung, (2) NIH3T3, (3) liver, (4) testis lysate; Predicted size: 54 kDa.



Western blot testing of Mkl antibody at 1:2000 dilution and mouse samples: (1) testis, (2) liver, (3) NIH3T3 lysate; Predicted size: 54 kDa.



Western blot testing of Mkl antibody at 1:1000 dilution and mouse samples: (1) lung, (2) spleen lysate; Predicted size: 54 kDa.

Description

The protein kinase domain is predicted to be catalytically inactive. Molecular function: protein binding. There are two isoforms.

Application Notes

Titration of the Mkl antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 444-472 from the mouse protein was used as the immunogen for this Mkl antibody.

Storage

Aliquot the Mkl antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.