


MLL5 Antibody (F40100)

Catalog No.	Formulation	Size
F40100-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40100-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8IZD2
Applications	Western Blot : 1:1000
Limitations	This MLL5 antibody is available for research use only.



Description

This gene is a member of the myeloid/lymphoid or mixed-lineage leukemia (MLL) family and encodes a protein with an N-terminal PHD zinc finger and a central SET domain. Overexpression of the protein inhibits cell cycle progression. Alternate transcriptional splice variants have been characterized. [provided by RefSeq].

Application Notes

Titration of the MLL5 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 93-120 from the human protein was used as the immunogen for this MLL5 antibody.

Storage

Aliquot the MLL5 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.