


MMP3 Antibody (R33577)


Catalog No.	Formulation	Size
R33577-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4314
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This MMP3 antibody is available for research use only.



Western blot testing of human prostate lysate with MMP3 antibody at 0.2ug/ml. Molecular weight: the expected ~54 kDa band and the additional ~24kDa band are both blocked by the immunizing peptide.



Western blot testing of mouse eye lysate with MMP3 antibody at 0.3ug/ml. Predicted molecular weight ~54 kDa.

Description

Additional name(s) for this target protein: Matrix metalloproteinase 3

Application Notes

Optimal dilution of the MMP3 antibody should be determined by the researcher.

Immunogen

Amino acids DPNAKKVHTLKSNS were used as the immunogen for this MMP3 antibody.

Storage

Aliquot and store the MMP3 antibody at -20oC.