

MMP7 Antibody (R31469)


Catalog No.	Formulation	Size
R31469	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with Trehalose
Gene ID	17393
Applications	Western Blot : 0.5-1ug/ml IHC (FFPE) : 0.5-1ug/ml ELISA : 0.1-0.5ug/ml (mouse protein tested); request BSA free format for coating
Limitations	This MMP7 antibody is available for research use only.



Western blot testing of MMP7 antibody and rat brain lysate. Predicted molecular weight ~30kDa.



IHC-P: MMP7 antibody testing of mouse kidney tissue



IHC-P: MMP7 antibody testing of rat intestine tissue

Description

Matrix metalloproteinase-7, previously called putative metalloproteinase I (PUMP1) or Matrilysin has been identified through studies of collagenase-related connective-tissue-degrading metalloproteinases produced by human tumors. The MMP7 protein has 267 amino acids and is significantly shorter than stromelysin or collagenase (477 and 469 amino acids, respectively). Matrix metalloproteinases play a crucial role in tumor invasion and metastasis. MMP7, a member of the matrix metalloproteinase family, is structurally different from the other matrix metalloproteinases by virtue of the absence of a conserved COOH-terminal protein domain. In addition, MMP7 mRNA is regulated in a specific and distinct manner in normal and malignant tissues. MMP7 has been shown to correlate with nodal or distant metastasis in colorectal carcinomas, however, its implication in early invasive colorectal carcinomas has not been determined. MMP7 is also a mediator of pulmonary fibrosis and a potential therapeutic target.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the MMP7 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Mouse partial recombinant protein (AA 95-264) was used as the immunogen for this MMP7 antibody.

Storage

After reconstitution, the MMP7 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.