

MMP7 Antibody (R35332)


Catalog No.	Formulation	Size
R35332-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4316
Applications	Western Blot : 0.03-0.1ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This MMP7 antibody is available for research use only.



Western blot testing of human breast lysate with MMP7 antibody at 0.03ug/ml. Predicted molecular weight ~30kDa.



IHC staining of FFPE human uterus with MMP7 antibody at 3.8ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Matrix metalloproteinase 7

Application Notes

Optimal dilution of the MMP7 antibody should be determined by the researcher.

Immunogen

Amino acids QKLYGKRSNSRKK were used as the immunogen for this MMP7 antibody.

Storage

Aliquot and store the MMP7 antibody at -20oC.