


## Moesin Antibody (R34765)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R34765-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                          |
|----------------------|----------------------------------------------------------|
| Availability         | 1-3 business days                                        |
| Species Reactivity   | Human, Mouse                                             |
| Predicted Reactivity | Rat, Dog                                                 |
| Format               | Antigen affinity purified                                |
| Clonality            | Polyclonal (goat origin)                                 |
| Isotype              | Goat Ig                                                  |
| Purity               | Antigen affinity                                         |
| Gene ID              | 4478                                                     |
| Applications         | Western Blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:16000 |
| Limitations          | This Moesin antibody is available for research use only. |



Western blot testing of human spleen lysate with Moesin antibody at 1ug/ml. Predicted molecular weight ~68 kDa but routinely observed at 68-78 kDa.

## Description

Additional name(s) for this target protein: MSN

## Application Notes

Optimal dilution of the Moesin antibody should be determined by the researcher.

## Immunogen

Amino acids NARDESKKTAND were used as the immunogen for this Moesin antibody.

## Storage

Aliquot and store the Moesin antibody at -20oC.