


MRC1L1 Antibody / MRC1 / Macrophage mannose receptor 1 (F40105)


Catalog No.	Formulation	Size
F40105-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40105-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P22897
Applications	Western Blot : 1:1000
Limitations	This MRC1L1 antibody is available for research use only.



Western blot testing of MRC1L1 antibody at 1:2000 dilution and human samples: Lane 1: lung lysate; 2: spleen lysate; 3: placenta lysate; Predicted band size : 166 kDa.



MRC1L1 antibody western blot analysis in mouse lung tissue lysate

Description

This gene encodes a type I membrane receptor protein highly related to the mannose receptor C type 1 protein (MRC1). The two genes are located near each other on chromosome 10 in a region thought to represent a segmental duplication.

Application Notes

Titration of the MRC1L1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 359-388 from the human protein was used as the immunogen for this MRC1L1 antibody.

Storage

Aliquot the MRC1L1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.