


MTHFD1 Antibody / C-1-tetrahydrofolate synthase (R35615)

Catalog No.	Formulation	Size
R35615-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4522
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This MTHFD1 antibody is available for research use only.



Western blot testing of human substantia nigra lysate with MTHFD1 antibody at 1ug/ml.
Predicted molecular weight: ~102 kDa.

Description

Additional name(s) for this target protein: Methylenetetrahydrofolate dehydrogenase (NADP+ dependent) 1

Application Notes

Optimal dilution of the MTHFD1 antibody should be determined by the researcher.

Immunogen

Amino acids RGDLNDCFIPCTPK were used as the immunogen for this MTHFD1 antibody.

Storage

Aliquot and store the MTHFD1 antibody at -20oC.