

MUC3B Antibody (F41684)


Catalog No.	Formulation	Size
F41684-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F41684-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9H195
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This MUC3B antibody is available for research use only.

250
130
95
72
55

MUC3B antibody western blot analysis in NCI-H292 lysate



MUC3B antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human stomach tissue.

Description

Major glycoprotein component of a variety of mucus gels. Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces (By similarity).

Application Notes

Titration of the MUC3B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 866-895 from the human protein was used as the immunogen for this MUC3B antibody.

Storage

Aliquot the MUC3B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.