


MUNC18-3 Antibody (R35945)

Catalog No.	Formulation	Size
R35945-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6814
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This MUNC18-3 antibody is available for research use only.



Western blot testing of mouse [A] and rat [B] spleen lysates with MUNC18-3 antibody at 1ug/ml. Predicted molecular weight: ~68kDa.

Description

Additional name(s) for this target protein: STXBP3, syntaxin binding protein 3, MUNC18C

Application Notes

Optimal dilution of the MUNC18-3 antibody should be determined by the researcher.

Immunogen

Amino acids SKSIRRCKEINIS were used as the immunogen for this MUNC18-3 antibody.

Storage

Aliquot and store the MUNC18-3 antibody at -20oC.