


MYB Antibody [clone 1279CT309.289] (F52454)

Catalog No.	Formulation	Size
F52454-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52454-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, k
Clone Name	1279CT309.289
Purity	Purified
UniProt	P10242
Applications	Western Blot : 1:1000
Limitations	This MYB antibody is available for research use only.



Western blot analysis of lysate from CEM cell line using MYB antibody at 1:1000.
Predicted molecular weight ~72 kDa.

Description

Transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Plays an important role in the control of proliferation and differentiation of hematopoietic progenitor cells.

Application Notes

Titration of the MYB antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This MYB antibody was produced from a mouse immunized with a KLH conjugated synthetic peptide between 1-185 amino acids from the Central region of human MYB.

Storage

Aliquot the MYB antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.