


MYBL2 Antibody (F53830)

Catalog No.	Formulation	Size
F53830-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P10244
Applications	Western Blot : 1:1000-2000
Limitations	This MYBL2 antibody is available for research use only.



Western blot testing of human 1) Raji and 2) HeLa cell lysate with MYBL2 antibody at 1:2000. Predicted molecular weight: ~79 kDa, may also be observed at 100~110 kDa.

Description

Transcription factor involved in the regulation of cell survival, proliferation, and differentiation. Transactivates the expression of the CLU gene.

Application Notes

Titration of the MYBL2 antibody may be required due to differences in protocols and substrate sensitivity.

Immunogen

A portion of amino acids 245-278 from the human protein was used as the immunogen for the MYBL2 antibody.

Storage

Aliquot the MYBL2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.