


MYRIP Antibody (R34377)


Catalog No.	Formulation	Size
R34377-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	25924
Localization	Cytoplasmic
Applications	Western Blot : 0.3-1ug/ml IHC (FFPE) : 3.75ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This MYRIP antibody is available for research use only.




Western blot testing of human amygdala lysate with MYRIP antibody at 0.3ug/ml.
Predicted molecular weight: ~95 kDa.



IHC testing of FFPE human cortex (brain) tissue with MYRIP antibody at 3.75ug/ml.
Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.



IHC testing of FFPE human colon tissue with MYRIP antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.



IHC testing of FFPE human breast tissue with MYRIP antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Myosin VIIA and Rab interacting protein, SLAC2C

Application Notes

Optimal dilution of the MYRIP antibody should be determined by the researcher.

Immunogen

Amino acids KDLMEPALESVMY were used as the immunogen for this MYRIP antibody.

Storage

Aliquot and store the MYRIP antibody at -20oC.