


## NAC1 Antibody / BTBD14B (R34663)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R34663-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                             |
|----------------------|-------------------------------------------------------------|
| Availability         | 1-3 business days                                           |
| Species Reactivity   | Human                                                       |
| Predicted Reactivity | Dog                                                         |
| Format               | Antigen affinity purified                                   |
| Clonality            | Polyclonal (goat origin)                                    |
| Isotype              | Goat Ig                                                     |
| Purity               | Antigen affinity                                            |
| Gene ID              | 112939                                                      |
| Applications         | Western Blot : 0.5-1.5ug/ml<br>ELISA (peptide) LOD : 1:2000 |
| Limitations          | This NAC1 antibody is available for research use only.      |



Western blot testing of human bone marrow lysate with NAC1 antibody at 0.5ug/ml.  
Predicted molecular weight: ~57 kDa, commonly observed at 57-70 kDa.

## Description

Additional name(s) for this target protein: BTB (POZ) domain containing 14B; BTBD14B

## Application Notes

Optimal dilution of the NAC1 antibody should be determined by the researcher.

## Immunogen

Amino acids KTEQQESDSVQC were used as the immunogen for this NAC1 antibody.

## Storage

Aliquot and store the NAC1 antibody at -20oC.