


## NAT16 Antibody / C7orf52 (F53914)

| Catalog No.   | Formulation                                | Size    |
|---------------|--------------------------------------------|---------|
| F53914-0.2ML  | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.2 ml  |
| F53914-0.05ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.05 ml |

[Bulk quote request](#)

|                           |                                                         |
|---------------------------|---------------------------------------------------------|
| <b>Availability</b>       | 1-3 business days                                       |
| <b>Species Reactivity</b> | Human                                                   |
| <b>Format</b>             | Antigen affinity purified                               |
| <b>Clonality</b>          | Polyclonal (rabbit origin)                              |
| <b>Isotype</b>            | Rabbit Ig                                               |
| <b>Purity</b>             | Antigen affinity                                        |
| <b>UniProt</b>            | Q8N8M0                                                  |
| <b>Applications</b>       | Western Blot : 1:1000-2000                              |
| <b>Limitations</b>        | This NAT16 antibody is available for research use only. |



Western blot testing of human A549 cell lysate with NAT16 antibody at 1:2000. Predicted molecular weight: 41 kDa.

### Description

Alternate names: n-Acetyltransferase 16

### Application Notes

Titration of the NAT16 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

A portion of amino acids 12-46 from the human protein was used as the immunogen for the NAT16 antibody.

### Storage

Aliquot the NAT16 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.