

NCAPH Antibody (R36151)

Catalog No.	Formulation	Size
R36151-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23397
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NCAPH antibody is available for research use only.



Western blot testing of K562 cell lysate with NCAPH antibody at 0.01ug/ml. Predicted molecular weight: ~83kDa, routinely observed at 90~97kDa.

Description

Additional name(s) for this target protein: Non-SMC condensin I complex, subunit H; HCAP-H

Application Notes

Optimal dilution of the NCAPH antibody should be determined by the researcher.

Immunogen

Amino acids TEDLSDVLVRQGD were used as the immunogen for this NCAPH antibody.

Storage

Aliquot and store the NCAPH antibody at -20oC.

References (1)