


NDUFS7 Antibody (R36039)

Catalog No.	Formulation	Size
R36039-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	374291
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This NDUFS7 antibody is available for research use only.



Western blot testing of human heart lysate with NDUFS7 antibody at 0.1ug/ml. Predicted molecular weight: ~24kDa but routinely observed at ~20kDa.

Description

Additional name(s) for this target protein: NADH dehydrogenase (ubiquinone) Fe-S protein 7

Application Notes

Optimal dilution of the NDUFS7 antibody should be determined by the researcher.

Immunogen

Amino acids SRGEYVAKLD were used as the immunogen for this NDUFS7 antibody.

Storage

Aliquot and store the NDUFS7 antibody at -20oC.