

NEDD4 Antibody (F52006)


Catalog No.	Formulation	Size
F52006-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52006-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	P46934
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:50-1:100 Flow Cytometry : 1:10-1:50
Limitations	This NEDD4 antibody is available for research use only.

250
130
95
72

Western blot analysis of NEDD4 antibody and mouse NIH3T3 lysate. Predicted molecular weight ~149 kDa; commonly observed at 110-149 kDa with a possible ~95 kDa cleavage band.



NEDD4 antibody IHC analysis in formalin fixed and paraffin embedded brain tissue.



NEDD4 antibody flow cytometric analysis of NCI-H460 cells (green) compared to a [negative control](http://search_result.php?search_txt=n1001) (blue). FITC-conjugated goat-anti-rabbit secondary Ab was used for the analysis.

Description

Nedd4 is one of a group of mouse genes that show developmentally regulated expression in mouse embryonic brain. Nedd4 is expressed in various other embryonic tissues and persists in most adult tissues.

Application Notes

Titration of the NEDD4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1288-1319 from the human protein was used as the immunogen for this NEDD4 antibody.

Storage

Aliquot the NEDD4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.