


Neurofascin Antibody (R34796)

Catalog No.	Formulation	Size
R34796-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23114
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Neurofascin antibody is available for research use only.



Western blot testing of human lung lysate with Neurofascin antibody at 0.5ug/ml.
Molecular weight: observed here at ~170kDa.

Description

Additional name(s) for this target protein: NFASC

Application Notes

Optimal dilution of the Neurofascin antibody should be determined by the researcher.

1. This Neurofascin antibody will detect isoforms 1/2/3/4.

Immunogen

Amino acids NFSPNQTKFTVQRT were used as the immunogen for this Neurofascin antibody.

Storage

Aliquot and store the Neurofascin antibody at -20oC.