

Neurofibromin 2 Antibody (Zebrafish) (R33848)

Catalog No.	Formulation	Size
R33848-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q32PS2
Applications	Immunohistochemistry (FFPE) : suitable (1) ELISA (peptide) LOD : 1:32000
Limitations	This Neurofibromin 2 antibody is available for research use only.



Description

Additional name(s) for this target protein: NF2, NF2a

Application Notes

Optimal dilution of the Neurofibromin 2 antibody should be determined by the researcher.

1. Reported to detect Neurofibromin 2 in FFPE zebrafish brain.

Immunogen

Amino acids EKSKHLQDQLK were used as the immunogen for this Neurofibromin 2 antibody.

Storage

Aliquot and store the Neurofibromin 2 antibody at -20oC.