


Neurofilament Antibody (light chain) [clone MSI27-1] (V7122)

Catalog No.	Formulation	Size
V7122-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V7122-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V7122SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug
V7122IHC-7ML	Prediluted in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide; *For IHC use only*	7 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	MSI27-1
Purity	Protein G affinity chromatography
UniProt	P07196
Localization	Cytoplasmic, membranous
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT Prediluted IHC Only Format : incubate for 30 min at RT (1)
Limitations	This Neurofilament antibody is available for research use only.



IHC testing of FFPE human cerebellum with Neurofilament antibody (clone MSI27-1). Staining of formalin-fixed tissues requires boiling tissue sections in pH 9 10mM Tris with 1mM EDTA for 10-20 min followed by cooling at RT for 20 min.

Description

Neurofilaments are type IV intermediate filament heteropolymers composed of light, medium, and heavy chains. Neurofilaments comprise the axoskeleton and they functionally maintain the neuronal caliber. They may also play a role in intracellular transport to axons and dendrites. [RefSeq]

Application Notes

Optimal dilution of the Neurofilament antibody should be determined by the researcher.

1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Immunogen

Recombinant human neurofilament light chain protein was used as the immunogen for the Neurofilament antibody.

Storage

Store the Neurofilament antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).