


Neurologin 3 Antibody (R35644)


Catalog No.	Formulation	Size
R35644-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54413
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Neurologin 3 antibody is available for research use only.



Western blot testing of human brain (amygdala) lysate with Neurologin 3 antibody at 0.1ug/ml. Predicted molecular weight: ~91kDa.



Western blot of HEK293 lysate overexpressing Neuregulin 3 probed with Neuregulin 3 antibody (mock transfection in lane 1). Predicted molecular weight: ~91kDa.

Description

Additional name(s) for this target protein: NLGN3

Application Notes

Optimal dilution of the Neuroligin 3 antibody should be determined by the researcher.

Immunogen

Amino acids RKDKRRQEPLRQPS were used as the immunogen for this Neuroligin 3 antibody.

Storage

Aliquot and store the Neuroligin 3 antibody at -20oC.