


Neuromedin-S Antibody / Nms (RQ7098)


Catalog No.	Formulation	Size
RQ7098	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q5H8A1
Applications	Immunohistochemistry (FFPE) : 2-5ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This Neuromedin-S antibody is available for research use only.



IHC staining of FFPE rat brain tissue with Neuromedin-S antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE mouse brain tissue with Neuromedin-S antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.

Description

Neuromedin S is a 36-amino acid neuropeptide found in the brain of humans and other mammals. This gene encodes a member of the neuromedin family of neuropeptides. The encoded protein is a precursor that is proteolytically processed to generate a biologically active neuropeptide that plays a role in the regulation of circadian rhythm, anorexigenic action, antidiuretic action, cardiovascular function and stimulation of oxytocin and vasopressin release. Mice lacking the encoded neuropeptide exhibit decreased heart rate without any accompanying changes in blood pressure. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar processing to generate the mature peptide.

Application Notes

Optimal dilution of the Neuromedin-S antibody should be determined by the researcher.

Immunogen

Recombinant mouse protein (amino acids D36-N144) was used as the immunogen for the Neuromedin-S antibody.

Storage

After reconstitution, the Neuromedin-S antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.