


NF-M Antibody [Neurofilament, medium] (R35201)

Catalog No.	Formulation	Size
R35201-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4741
Localization	Cytoplasmic, membranous
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NF-M antibody is available for research use only.



Western blot testing of mouse brain lysate with NF-M antibody at 0.03ug/ml. Predicted molecular weight: ~102kDa but routinely observed at 145~160kDa.

Description

Additional name(s) for this target protein: Neurofilament medium, NEFM

Application Notes

Optimal dilution of the NF-M antibody should be determined by the researcher.

Immunogen

Amino acids QASHITVERKDYLK were used as the immunogen for this NF-M antibody.

Storage

Aliquot and store the NF-M antibody at -20oC.