


NFAT2 Antibody / NFATC1 (isoforms A/B/C) (R32829)


Catalog No.	Formulation	Size
R32829	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	O95644
Localization	Cytoplasmic, nuclear
Applications	Western Blot : 0.5-1ug/ml Flow Cytometry : 1-3ug/million cells Direct ELISA (human Recombinant Protein : 0.1-0.5ug/ml (BSA-free formulation available)
Limitations	This NFAT2 antibody is available for research use only.



Flow cytometry testing of fixed and permeabilized human Jurkat cells with NFAT2 antibody at 1ug/10⁶ cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue=NFAT2 antibody.



Western blot testing of human 1) Jurkat, 2) Daudi, 3) Ramos and 4) Raji cell lysate with NFAT2 antibody at 0.5ug/ml. Predicted molecular weight: ~77 kDa (isoform A), ~88 kDa (isoform B), ~101 kDa (isoform C).

