


NFKBIE Antibody (F53341)

Catalog No.	Formulation	Size
F53341-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53341-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O00221
Applications	Western Blot : 1:2000
Limitations	This NFKBIE antibody is available for research use only.



Western blot testing of NFKBIE antibody at 1:2000 dilution. Lane 1: PC3 lysate; 2: A375 lysate; 3: HeLa lysate; 4: LNCaP lysate; 5: THP-1 lysate; Predicted band size : 53 kDa.

Description

Inhibits NF-kappa-B by complexing with and trapping it in the cytoplasm. Inhibits DNA-binding of NF-kappa-B p50-p65 and p50-c-Rel complexes.

Application Notes

Titration of the NFKBIE antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This NFKBIE antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 109-141 amino acids from the N-terminal region of human NFKBIE.

Storage

Aliquot the NFKBIE antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.