


NGAL Antibody (R34280)

Catalog No.	Formulation	Size
R34280-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	16819
Localization	Cytoplasmic
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This NGAL antibody is available for research use only.



Western blot testing of mouse kidney lysate with NGAL antibody at 2ug/ml. Predicted molecular weight: 22-25kDa.

Description

Additional name(s) for this target protein: Lipocalin 2; Lcn2

Application Notes

Optimal dilution of the NGAL antibody should be determined by the researcher.

Immunogen

Amino acids HRYPQVQSYNVQ were used as the immunogen for this NGAL antibody.

Storage

Aliquot and store the NGAL antibody at -20oC.