


NGFR Antibody (R34988)

Catalog No.	Formulation	Size
R34988-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4804
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This NGFR antibody is available for research use only.



Western blot testing of rat heart lysate with NGFR antibody at 1ug/ml. Predicted molecular weight ~45/75kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Nerve growth factor receptor, p75NTR, CD271

Application Notes

Optimal dilution of the NGFR antibody should be determined by the researcher.

Immunogen

Amino acids PEGSDSTAPSTQEP were used as the immunogen for this NGFR antibody. The amino acid sequence used

as immunogen is from the extracellular portion of the protein.

Storage

Aliquot and store the NGFR antibody at -20oC.