

NKX1.2 Antibody (F41662)


Catalog No.	Formulation	Size
F41662-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F41662-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9UD57
Applications	Western Blot : 1:1000 Flow Cytometry : 1:10-1:50
Limitations	This NKX1.2 antibody is available for research use only.



NKX1.2 antibody western blot analysis in 293 lysate



NKX1.2 antibody flow cytometric analysis of 293 cells (green) compared to a [negative control](#) (blue). FITC-conjugated goat-anti-rabbit secondary Ab was used for the analysis.

Description

NKX1-2 may function in cell specification, particularly in the CNS (By similarity).

Application Notes

Titration of the NKX1.2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-2 from the human protein was used as the immunogen for this NKX1.2 antibody.

Storage

Aliquot the NKX1.2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.