


NOS1AP Antibody (R35153)

Catalog No.	Formulation	Size
R35153-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9722
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This NOS1AP antibody is available for research use only.



IHC testing of FFPE human brain with NOS1AP antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: Nitric oxide synthase 1 adaptor protein, CAPON

Application Notes

Optimal dilution of the NOS1AP antibody should be determined by the researcher.

Immunogen

Amino acids PSKTKYNLVDDGH were used as the immunogen for this NOS1AP antibody.

Storage

Aliquot and store the NOS1AP antibody at -20oC.