


NOTCH3 Antibody (R35352)

Catalog No.	Formulation	Size
R35352-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4854
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This NOTCH3 antibody is available for research use only.



IHC testing of FFPE human hippocampus with NOTCH3 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: CADASIL, CASIL

Application Notes

Optimal dilution of the NOTCH3 antibody should be determined by the researcher.

Immunogen

Amino acids QLGPQPEVTPKRQ were used as the immunogen for this NOTCH3 antibody.

Storage

Aliquot and store the NOTCH3 antibody at -20oC.