


Npc2 Antibody / Niemann pick c2 (R35185)

Catalog No.	Formulation	Size
R35185-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	67963
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Npc2 antibody is available for research use only.



Western blot testing of mouse kidney lysate with Npc2 antibody at 0.1ug/ml.
Predicted/observed molecular weight: ~17 kDa.

Description

Additional name(s) for this target protein: Niemann-Pick disease, type C2

Application Notes

Optimal dilution of the Npc2 antibody should be determined by the researcher.

Immunogen

Amino acids QAEPLHFKDCGSK were used as the immunogen for this Npc2 antibody.

Storage

Aliquot and store the Npc2 antibody at -20oC.