


NPC2 Antibody / Niemann Pick C2 (R35186)

Catalog No.	Formulation	Size
R35186-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10577
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This NPC2 antibody is available for research use only.



Western blot testing of human testis lysate with NPC2 antibody at 0.5ug/ml. Predicted molecular weight: ~17 kDa, can be observed as a ~21/23 kDa doublet in some samples. (Ref 1).

Description

Additional name(s) for this target protein: Niemann-Pick disease, type C2

Application Notes

Optimal dilution of the NPC2 antibody should be determined by the researcher.

Immunogen

Amino acids QAEPVQFKDCGSVD were used as the immunogen for this NPC2 antibody.

Storage

Aliquot and store the NPC2 antibody at -20oC.

References (1)