


NPFFR1 Antibody (R35720)

Catalog No.	Formulation	Size
R35720-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64106
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This NPFFR1 antibody is available for research use only.



Western blot testing of human cerebellum lysate with NPFFR1 antibody at 0.03ug/ml.
Predicted/observed molecular weight: ~48kDa.

Description

Additional name(s) for this target protein: Neuropeptide FF receptor 1

Application Notes

Optimal dilution of the NPFFR1 antibody should be determined by the researcher.

Immunogen

Amino acids RPSGSHKEAYSERP were used as the immunogen for this NPFFR1 antibody.

Storage

Aliquot and store the NPFFR1 antibody at -20oC.