


## Npy2r Antibody (R34349)

Catalog No.	Formulation	Size
R34349-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4887
Applications	Western Blot : 0.1-0.3ug/ml Immunoprecipitation : suitable ELISA (peptide) LOD : 1:32000
Limitations	This Npy2r antibody is available for research use only.



Western blot testing of rat brain lysate with Npy2r antibody at 0.1ug/ml. Predicted molecular weight: ~43 kDa.

## Description

Additional name(s) for this target protein: Neuropeptide Y receptor Y2

## Application Notes

Optimal dilution of the Npy2r antibody should be determined by the researcher.

## Immunogen

Amino acids KAKKNLEVKKNN were used as the immunogen for this Npy2r antibody.

## **Storage**

Aliquot and store the Npy2r antibody at -20oC.