


## NR0B1 Antibody / DAX-1 (R33235)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33235-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                              |
|--------------------|--------------------------------------------------------------|
| Availability       | 1-3 business days                                            |
| Species Reactivity | Human                                                        |
| Format             | Antigen affinity purified                                    |
| Clonality          | Polyclonal (goat origin)                                     |
| Isotype            | Goat Ig                                                      |
| Purity             | Antigen affinity                                             |
| Gene ID            | 190                                                          |
| Applications       | Western Blot : 0.5-1.5ug/ml<br>ELISA (peptide) LOD : 1:32000 |
| Limitations        | This NR0B1 antibody is available for research use only.      |



Western blot testing of human testis lysate with NR0B1 antibody at 0.5ug/ml. Predicted molecular weight: 44-52 kDa (two isoforms).

### Description

Additional name(s) for this target protein: Nuclear receptor subfamily 0, group B, member 1; DAX1; Dosage-sensitive sex reversal, adrenal hypoplasia critical region, on chromosome X, gene 1

### Application Notes

Optimal dilution of the NR0B1 antibody should be determined by the researcher.

### Immunogen

Amino acids CGEDHPQQGSTLY were used as the immunogen for this NR0B1 antibody.

## Storage

Aliquot and store the NR0B1 antibody at -20oC.