


NR2F1 Antibody / COUP TF1 [clone ACHO-14] (RQ5326)

Catalog No.	Formulation	Size
RQ5326	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	ACHO-14
Purity	Affinity purified
UniProt	P10589
Applications	Western Blot : 1:500-1:1000
Limitations	This NR2F1 antibody is available for research use only.



Western blot testing of human HEK293 cell lysate with NR2F1 antibody. Predicted molecular weight ~46 kDa.

Description

Coup (chicken ovalbumin upstream promoter) transcription factor, also called Nuclear receptor subfamily 2 group F member 1 or NR2F1, binds to the ovalbumin promoter and, in conjunction with another protein (S300-II) stimulates initiation of transcription. Binds to both direct repeats and palindromes of the 5'-AGGTCA-3' motif. Represses transcriptional activity of LHCG. [UniProt]

Application Notes

Optimal dilution of the NR2F1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human COUP TF1 / NR2F1 was used as the immunogen for the NR2F1 antibody.

Storage

Store the NR2F1 antibody at -20oC.