


NR5A2 Antibody (R36406)

Catalog No.	Formulation	Size
R36406-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2494
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NR5A2 antibody is available for research use only.



Western blot testing of HepG2 lysate with NR5A2 antibody at 1ug/ml. Predicted molecular weight: ~56/61kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: Liver receptor homolog 1 (LRH1); nuclear receptor subfamily 5, group A, member 2

Application Notes

Optimal dilution of the NR5A2 antibody should be determined by the researcher.

1. This sequence is common to isoforms 1 and 2.

Immunogen

Amino acids QEQANRSKHEK were used as the immunogen for this NR5A2 antibody. (1)

Storage

Aliquot and store the NR5A2 antibody at -20oC.