


## NRAS Antibody [clone 822CT17.5.1] (F52385)

Catalog No.	Formulation	Size
F52385-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52385-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Chicken, Mouse, Pig
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1
Clone Name	822CT17.5.1
Purity	Purified
UniProt	P01111
Applications	Western Blot : 1:1000
Limitations	This NRAS antibody is available for research use only.



NRAS antibody western blot analysis in 1) human A431 and 2) rat C6 lysate. Predicted molecular weight ~21 kDa.

### Description

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

### Application Notes

Titration of the NRAS antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

A portion of amino acids 70-101 from the human protein was used as the immunogen for this NRAS antibody.

## **Storage**

Aliquot the NRAS antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.