


NRAS Antibody (F54081)

Catalog No.	Formulation	Size
F54081-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54081-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Rat
Predicted Reactivity	Human, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P01111
Applications	Western Blot : 1:2000
Limitations	This NRAS antibody is available for research use only.



Western blot testing of NRAS antibody at 1:2000 + rat PC12 cell lysate. Predicted molecular weight ~21 kDa.

Description

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the NRAS antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Recombinant human NRAS was used as the immunogen for the NRAS antibody.

Storage

Aliquot the NRAS antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.