


NUDC Antibody (R35084)

Catalog No.	Formulation	Size
R35084-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10726
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This NUDC antibody is available for research use only.



Western blot testing of HeLa lysate with NUDC antibody at 0.3ug/ml. Predicted molecular weight: ~38kDa.

Description

Additional name(s) for this target protein: Nuclear distribution gene C

Application Notes

Optimal dilution of the NUDC antibody should be determined by the researcher.

Immunogen

Amino acids PNLGNGADLPNYR were used as the immunogen for this NUDC antibody.

Storage

Aliquot and store the NUDC antibody at -20oC.