


OCT1 Antibody / POU2F1 (R34943)

Catalog No.	Formulation	Size
R34943-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5451
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This OCT1 antibody is available for research use only.



Western blot testing of nuclear HeLa lysate with OCT1 antibody at 0.5ug/ml. Predicted molecular weight: ~77 kDa, routinely observed at 90-100 kDa.

Description

Additional name(s) for this target protein: Octamer-binding protein 1; Octamer-binding transcription factor 1; POU domain, class 2, transcription factor 1

Application Notes

Optimal dilution of the OCT1 antibody should be determined by the researcher.

Immunogen

Amino acids NVQSKSNEESGDSQQ were used as the immunogen for this OCT1 antibody.

Storage

Aliquot and store the OCT1 antibody at -20oC.