


OR2AE1 Antibody / C-Terminal Region (F54216)

Catalog No.	Formulation	Size
F54216-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54216-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal
Isotype	Rabbit Ig
Purity	Peptide affinity purified
UniProt	Q8NHA4
Applications	Western Blot : 1:1000
Limitations	This OR2AE1 antibody is available for research use only.



Western blot testing of human 1) U-87 MG, 2) U-2 OS and 3) MDA-MB-231 cell lysate with OR2AE1 antibody at 1:1000. Predicted molecular weight ~37 kDa.

Description

Olfactory receptor 2AE1 is an odorant receptor.

Application Notes

The stated application concentrations are suggested starting points. Titration of the OR2AE1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 278-312 from the C-terminal region of the human protein was used as the immunogen for the OR2AE1 antibody.

Storage

Aliquot the OR2AE1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.