


OR5V1 Antibody / C-Terminal Region (F54212)

Catalog No.	Formulation	Size
F54212-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54212-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal
Isotype	Rabbit Ig
Purity	Peptide affinity purified
UniProt	Q9UGF6
Applications	Western Blot : 1:1000
Limitations	This OR5V1 antibody is available for research use only.



Western blot testing of human 1) SH-SY5Y, 2) PANC-1 and 3) MCF-7 cell lysate with OR5V1 antibody at 1:1000. Predicted molecular weight ~36 kDa.

Description

Olfactory receptor 5V1 is an odorant receptor.

Application Notes

The stated application concentrations are suggested starting points. Titration of the OR5V1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 241-274 from the C-terminal region of the human protein was used as the immunogen for the OR5V1 antibody.

Storage

Aliquot the OR5V1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.