


OSBPL9 Antibody / ORP9 (R33918)

Catalog No.	Formulation	Size
R33918-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	114883
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This OSBPL9 antibody is available for research use only.



IHC testing of FFPE human lung with OSBPL9 antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Oxysterol binding protein-like 9; OSBP-related protein 9; ORP9

Application Notes

Optimal dilution of the OSBPL9 antibody should be determined by the researcher.

Immunogen

Amino acids EPLLKRLGAAKH were used as the immunogen for this OSBPL9 antibody.

Storage

Aliquot and store the OSBPL9 antibody at -20oC.