


OSMR Antibody (R35303)

Catalog No.	Formulation	Size
R35303-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9180
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This OSMR antibody is available for research use only.



Western blot testing of HeLa lysate with OSMR antibody at 0.1ug/ml. Predicted molecular weight: ~111kDa.

Description

Additional name(s) for this target protein: Oncostatin M receptor

Application Notes

Optimal dilution of the OSMR antibody should be determined by the researcher.

Immunogen

Amino acids QIHGEQLDPHVT were used as the immunogen for this OSMR antibody. The amino acid sequence used as immunogen is from the extracellular portion of the protein.

Storage

Aliquot and store the OSMR antibody at -20oC.