


## OXER1 Antibody / Oxoeicosanoid receptor 1 (R36055)

Catalog No.	Formulation	Size
R36055-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	165140
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This OXER1 antibody is available for research use only.



Western blot testing of human lung lysate with OXER1 antibody at 1ug/ml. Predicted molecular weight ~41 kDa.

### Description

Additional name(s) for this target protein: Oxoeicosanoid (OXE) receptor 1; 5-oxo-ETE GPCR

### Application Notes

Optimal dilution of the OXER1 antibody should be determined by the researcher.

### Immunogen

Amino acids SRQWRYREASRK were used as the immunogen for this OXER1 antibody.

## Storage

Aliquot and store the OXER1 antibody at -20oC.