


p16INK4a Antibody (R35748)


Catalog No.	Formulation	Size
R35748-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	12578
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This p16INK4a antibody is available for research use only.



Western blot testing of mouse liver lysate with p16INK4a antibody at 0.3ug/ml. Predicted molecular weight ~16 kDa.



Western blot testing of 1) mouse and 2) rat liver lysate with p16INK4a antibody at 0.3ug/ml. Predicted molecular weight ~16 kDa.

Description

Additional name(s) for this target protein: Cyclin-dependent kinase inhibitor 2A; CDKN2A; p19ARF; p16

Application Notes

Optimal dilution of the p16INK4a antibody should be determined by the researcher.

Immunogen

Amino acids RRGPHRNPGPGDDD were used as the immunogen for this p16INK4a antibody.

Storage

Aliquot and store the p16INK4a antibody at -20oC.